

Step-By-Step Tutorial: Enabling the Athletic Trainer System (ATS) Integration Support

If you need help enabling your integration please reach out to: support@xlntbrain.raiseaticket.com.

1. Upon logging into your XLNTBrain account, select Integrations on the left navigation
2. Athletic Trainer System will appear to have a status of “Not Configured” when not set up, to enable select Configure Athletic Trainer System

Integrations *Integration management for Xavier School for the Gifted*

Note: This page allows you to configure credentials and settings for different integrations with external systems. Make sure you have the correct credentials before making any changes.

Available Integrations

Healthy Roster

Integration for sending concussion reports to the Healthy Roster platform. This integration allows secure transmission of PDF reports and patient demographic data.

Status: Online

[Configure Healthy Roster](#)

Athletic Trainer System

Integration to send basic concussion data to Athletic Trainer System. This integration allows to upload Athletes Concussion data into ATS.

Status: Online

[Configure Athletic Trainer System](#)

More Integrations

More integration options with other systems will be added soon. Contact support to request new integrations.

Status: Coming Soon

[View options](#)

3. Review the instructions

Athletic Trainer System *Credentials for Xavier School for the Gifted*

[← Back to Integrations](#)

Instructions

This integration allows uploading concussion and athlete data to Athletic Trainer System.

Requirements:

- You must enter the ATS Database Name to be used. It should be the same text you see in the login screen in ATS, it's a text field with the label "Database".
- ATS will import the uploaded data every noon and you should be able to search the concussion and athlete data in ATS app after the import is finished.

Important: This information is confidential and should not be shared with unauthorized staff. Make sure to enter the correct values.

4. Input your ATS Database Name where you want the data to be uploaded

API Configuration

ATS Database Name:

Database Name field that can be seen in the ATS login form.

5. Check if desired to apply automatic synchronization with Athletic Trainer System


Auto Upload: Enable automatic synchronization with Athletic Trainer System
By enabling this option, all your organization concussion and athlete data will be automatically uploaded nightly to the ATS System.

6. Check to agree to the terms of the integration

Data Consent: By enabling this integration, you agree to transfer concussion and athlete data on the XLNTBrain platform to Athletic Trainer System. If you wish to stop this from happening, please disable this integration. If you have any questions/comments, please contact us at support@xlntbrain.raiseaticket.com.

7. Select Save Configuration

Save Configuration

8. To test if your integration is active, send a test file. You will find an  notification if the connection is ok.

Test Connection

You can use this form to test concussion and athlete data upload to a test SFTP Server.

